

Bench Test Report: Cell Culture Flask

Cell Culture



293A – Human Embryonic Kidney Epithelial

www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:
Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or
Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL Cell Culture Flask – 293A



- Human kidney epithelial (293A) cells cultured on SPL Cell culture flask for 72 hours with seeding density of 6×10^5 cells/flask.
- Cells were maintained in 37°C humidified condition, in 10 ml of DMEM, supplemented with 10% FBS and 1% antibiotics.
- Cells were tested negative for the presence of microbial.

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd.