

Sample Tube (#50221)

1. Description

SPL Sample Tubes are designed for storing various clinical samples. SPL Life Sciences provides durable, easy-to-use bottles to meet various demands of customers.

2. Features

- ❖ Tight seal with screw cap
- ❖ External graduations
- ❖ Self-standing type conical bottom
- ❖ Autoclavable

3. Product Range

| Cat. No. | Material (Body/Cap) | External Dimension (mm) | Total Volume (ml) | Sterile | Packaging |
|----------|---------------------|-------------------------|-------------------|---------|-----------|
| 50221 | PP / PE | 25.00 x 73.00 | 20.00 | - | 25 / 500 |

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs




7. Specifications & Technical Drawing

- Specifications

| Dimensions | Material (Resin) | Packaging |
|--|---|---|
| Diameter: 25.20 mm Height: 75.00 mm | ▪ Polypropylene (PP, body) ▪ High Density Polyethylene (HDPE, Cap) | ▪ Pieces/Sleeve: 25 ▪ Pieces/Case: 500 |

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:
Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr
To place an order, contact your local distributor or
Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-SAT

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 2, July 13th, 2016

