

Snap Tube Rack

(#53015/53050)

1. Description

SPL Snap Tube Racks have identical design to the SPL Conical Tube Racks, but has significant differences in inner dimensions and enhanced structural integrity. Racks are uniquely designed for handling and storage of 25 Snap Tubes. (5ml/15ml or 50ml respectively)

2. Features

- ❖ Excellent stackability
- ❖ Numeric labeling
- ❖ 2 colors : Blue, Green


3. Product Range

Cat. No.	Material	Description	Packaging
53015	PP	25 Holes for 5ml, 15 ml tubes	50
53050	PP	25 Holes for 50 ml tubes	20

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs



SPL-TDS-SNTURA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 2, July 25th, 2016


6. Specifications & Technical Drawing

- Specifications


Dimensions		Material (Resin)	Packaging	
53015	53050		53015	53050
Length: 113 mm	Length: 139 mm	Polypropylene (PP)	▪ Pieces/Case: 50	▪ Pieces/Case: 20
Width: 139.2 mm	Width: 262 mm			
Height: 48.1 mm	Height: 48.1 mm			

- Technical Drawing

Cat. No. 53015



Cat. No. 53050



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:
Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr
To place an order, contact your local distributor or
Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-SNTURA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 2, July 25th, 2016

