

PCR Tube Rack

(#80096)

1. Description

SPL PCR Tube Racks are especially designed to handle and store PCR Tubes or Strips.

2. Features

- ❖ 5 colors: White, Blue, Green, Orange, Purple.
- ❖ Autoclavable

3. Product Range

| Cat. No. | Material | Description | Packaging |
|----------|----------|-------------------|-----------|
| 80096 | PP | 8 x 12 (96 Holes) | 10 |

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs




6. Specifications & Technical Drawing

- Specifications

| Dimensions | | Material (Resin) | Packaging |
|------------|-----------|--------------------|-------------------|
| Length: | 86.00 mm | Polypropylene (PP) | ▪ Pieces/Case: 10 |
| Width: | 123.20 mm | | |
| Height: | 26.00 mm | | |

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-PCTURA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

© SPL Life Sciences Co., Ltd., Version 2, November 11th, 2016