

Biohazard Bag

(#98000/98001/98002)

1. Description

Biohazard Bag is designed for the safe storage, transportation and disposal of biohazardous waste. The Biohazard Bag is made of High Density Polyethylene and autoclave is possible, which ensure safely dispose of waste.

2. Features

- ❖ Thickness: 0.05 mm
- ❖ 3 different sizes
- ❖ Autoclavable
- ❖ Disposable

3. Product Range

Cat. No.	Material	Dimension (cm)	Packaging
98000	HDPE	30.00 x 45.00	100
98001	HDPE	40.00 x 60.00	100
98002	HDPE	60.00 x 90.00	100

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-BHB

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 1, December 21st, 2021